


3 Modes of Action:


Deepen Your Knowledge of the ToeFX Process



Mechanical:

The Debridement

Effectively and efficiently debride the toenail. Take time to break down the fungus build up on the nail. Keratin is hard to penetrate without excellent debridement.



Adhesive:

The Serum

The serum used in the ToeFX treatment is specially formulated to penetrate the entire nail, to promote adhesion to the fungal organisms that give rise to toenail fungus.



Physical:

The ClearToe

Therapy Light

The formula (penetration) and light (photodisinfection) work together to destroy the fungal cells.